



Lighting Cleaning Guide

Proper cleaning starts with the right products. Use only recommended cleaning and disinfecting solutions to help protect your Skytron lights, preserve their appearance, and ensure reliable performance over time.



Recommended Cleaning & Disinfection Guidelines

The ingredients in the following products have been reviewed and are considered compatible with Skytron lighting fixtures:

- Sani-Cloth AF3 Wipes
- Diversey Oxivir 1 Wipes
- Clorox Hydrogen Peroxide Wipes (1.4% hydrogen peroxide)
- Sani-Cloth HP1 Wipes (0.58% hydrogen peroxide)
- PREempt RTU
- PREempt Plus
- Ecolab Daily Disinfectant Cleaner

Skytron recommends consulting the disinfectant or cleaning solution manufacturer to verify compatibility with the following materials. To help prevent damage, avoid using abrasive cloths, paper towels, or scrubbing pads.

- Polycarbonate (PC)
- Polysulfone (PSU)
- Acrylic (PMMA)
- Silicone

After the required disinfectant contact time, wipe the diffuser with a soft, non-abrasive cloth to remove residue and prevent staining. When using UV disinfection, position the light toward the ceiling to avoid diffuser damage. Always follow the cleaner manufacturer's compatibility guidelines and avoid abrasive materials or cleaners containing high concentrations of alcohol, ethylene glycol, phenol, iodophors, or glutaraldehyde, as they may damage the fixture.

