



Lucina 4

Specialty Light



Keep Your Focus on the Patient


The Lucina 4 is perfect for bringing simple operation and flexible positioning to any clinical environment. This specialty light improves healthcare spaces by keeping them open and free of ceiling mounted equipment conflicts.

Ideal for:

- Diagnostic / Exam Light
- Labor, Delivery, Recovery and Postpartum Suites
- Intensive Care Units
- Hybrid Operating Rooms
- Cath Labs

Simple and Flexible Positioning

- Skytron exclusive strobe-guided technology
- Guide the light beam to the targeted location with a click of a button



Lucina 4 is a diagnostic/examination light and is not intended to be used as a primary surgical luminaire. Skytron recommends a surgical light (fixed or mobile) be included in clinical settings where surgical illumination may be required.



Simple Positioning



Clutter Free Environment



Cost Effective

Lucina 4 Details:

- Unique LED strobe-guided positioning wand
- Angles up to 47° for superior illumination
- Long lasting LEDs
- Fully recessed 24" square fixture
- Wall control with On/Off and intensity touchpad and (back-up) positioning buttons
- Available in one, two, three and four fixture configurations



Tranquil and Accessible Environment



Cost Effective

- Long-lasting life span with energy efficient LEDs
- Flush mounting within the ceiling, an ideal choice for low ceiling heights
- Clears the area and saves space for other equipment

Clean and Clear Environment

- Flush mounting structure to ensure that each clinical space remains comfortable and less intimidating for the patient
- Less invasive than traditional exam lighting arm configurations

LED Strobe-Guided Positioning Wand and Wall Mounted Intensity Control



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV4 03-24
PN L2-010-118