



Sterile Processing Department

Efficiency Solutions





Enhancing Sterile Processing Through Efficiency and Innovation

Skytron sterile processing equipment is engineered to optimize space, minimize water and chemical usage, and support the smooth, cost-efficient operation of your Sterile Processing Department. Our solutions are designed to boost productivity, reduce expenses, and ensure reliable, effective processing of medical devices.

Lower Cost of Ownership

- High quality and long lasting
- Rapid return on investment

Works Better with Less

- Use up to 63% less water
- Less chemicals
- Reduced electrical usage

Optimize Your Space

- Variety of sizes to choose from
- Process up to 37% more containers per cycle
- Work in noise levels as low as <63 dB



Water
Conservation



Lower Cost
of Ownership



Smaller
Footprint



Processing Sinks

Skytron's Processing Sinks compliment your first step in the disinfection process and provide a means for the pre-cleaning, inspection and sorting of surgical instruments.

Several Sizes Available

- Single, double, or triple basin options
- Universal workflow

Stationary or Electronic Height Adjustable Models

- 37" Standard height for stationary sinks
- Ergonomically friendly electronic adjustable height with four position programming
- Height adjustable settings ranging from 35"-42"

Clean Work Environment

- Laser etched gallon markings - quick and easy measurements
- 12" Backsplash to contain material spillage and optional backsplash mounted utilities
- Optional pegboard back for organized storage of cleaning accessories
- Included workspace for staging and inspection



Clean
Workspace



Comfortable Work
Environment



Flex Ultrasonic



Flex CT6 Countertop Ultrasonic

Ultrasonic Washer-Disinfectors

The Skytron Flex was designed with versatility and throughput in mind. This ultrasonic washer-disinfector has the ability to process large loads of non-proprietary trays, run mixed loads of da Vinci, orthopedic, cannulated, and non-cannulated instruments - helping you do more in less time! With Sterile Processing Departments facing increasing pressure to reduce costs and improve efficiency, the CT6 assists in achieving these goals by providing a reliable, effective, and affordable way to clean surgical instruments.

Flex Ultrasonic Features

- Triple the capacity of standard ultrasonic washer-disinfectors
- Two thermal disinfection cycles - Cannulated and Non-Cannulated
- Reprocess mixed loads in the same cycle
- Multi-step cycle cleans and disinfects in less than 35 minutes
- Meets AAMI Standards
- Utilize your current baskets and vendor trays
- With its thermal disinfection rinse, Flex can be used as a backup washer-disinfector

Flex CT6 Countertop Ultrasonic Features

- Compact, efficient, and eco-friendly
- Affordable ultrasonic solution
- Patented suction and ultrasonic cleaning technique
- Small environmental footprint, perfect for ASCs
- Exceeds the AAMI TIR 30 standard for cleaning effectiveness
- Consumes fewer utilities which contributes to environmental sustainability and cost savings



Clean Water &
Detergent Flush



Save Time
and Money



Reduce Labor
& Errors



Washer-Disinfectors

The Sterling M1000, P1000, P1800 and C2000 make up the Skytron family of Washer-Disinfectors. The Sterling portfolio provides a full array of throughput options to outfit any Sterile Processing Department, offering the lowest possible environmental impact. With a large capacity, small footprint, fast cycles, low water, detergent and electrical consumption, the Sterling series automates the decontamination process of surgical instruments.

Water Conservation

- Up to 65% less water than other washers
- Save 50% in chemical and usage costs

Comply With Your Medical Device IFUs

- Configurable cycle parameters to meet your needs
- Remote connectivity with cycle data storage

Boost Your Efficiency

- Best throughput per cycle in a 40% smaller footprint
- Recyclable water for more sustainable consumption



Save Water



Compliant



Smaller
Footprint



Cart Washer-Disinfector


The Sterling C2000 Cart Washer streamlines and automates the decontamination process — bringing unsurpassed efficiency and flexibility to the Sterile Processing Department. The C2000 is engineered to be environmentally friendly by conserving water, energy, and chemical consumption using its built-in control functions and water recirculation system.

Efficient Water Coverage Throughout the Chamber

- 20 wash arms with 120 oscillating nozzles
- Pressure controlled based on load configuration
- Wide spray pattern for complete coverage
- Spray arms and nozzles can be cleaned without tools
- Fits up to 3 carts in one washer
- Validated to disinfect instruments and loaners



Water, Chemical,
& Electricity
Conservation



Low
Noise
Level



Fast
Process
Time



ErgoStat® Prep & Pack Table



Stainless Prep & Pack Table

Prep & Pack Tables

The ErgoStat® Prep & Pack table and the Stainless Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly and packaging in your Sterile Processing Department. Consider your staff's well-being and choose a Prep & Pack table that easily adjusts to the various physical attributes and requirements of the human body – ultimately achieving better health, safety, and employee longevity.

ErgoStat® Features

- Increased stability with heavy duty inner/ outer leg extrusions
- 750 lb dynamic and 1,000 lb static weight capacity
- Skytron exclusive non-porous solid Corian or Stainless Surface available

Alleviate Physical Fatigue

- ErgoStat® - 26 ½" to 41 ½" height adjustable workspace
- ErgoStat® - 43" vertical riser frame raises and lowers with the base
- Ergonomically correct alignment of computer monitors, lighting, bins, and pegboards
- Powered and manual (pin-adjustable) models available



Built to Last



Alleviate
Physical Stress



Sterilizers

When it comes to sterilization technology, we understand that your needs are unique – in capacity, efficiency, utility usage savings, and more. Whether you're outfitting a new facility, renovating an older one, or simply replacing outdated technology, Integrity brings superior technology and versatility that improves your sterile processing workflow. We offer three different sizes of Integrity Sterilizers so we can provide your Sterile Processing Department with exactly what you need.

Work Efficiently in Small Spaces

- Integrity 215 has wider and taller chambers to accommodate larger trays horizontally
- Simple Plug-and-Play installation
- Save up to 63% in water consumption

Configurable and Easy to Use

- Electric power door offers hands-free control
- 12 preset standard cycles - 12 additional customizable cycles

Capacities to Meet Diverse Needs

- Integrity 175 models are compact - a great choice for a small department footprint
- Integrity 215 doubles competitive validated capacity from three 25lbs trays to six - that's 150lbs per cycle
- Integrity 270 processes up to 16 trays or 400lbs per cycle



Water
Conservation



Customizable
Cycles



Conserve
Space



EASE Modular Wall

Skytron's EASE modular walls offer a hygienic and healing design for any clinical environment. Our leading solution is great for surgical units but also suitable to all hospital departments that require easy to clean surfaces with flawless quality. Equipped with a self-supporting substructure and a variety of materials and finishes, the EASE wall is certain to set your facility apart from the rest.

Modular Wall Panels

- Resistant to UV, detergents, bactericidal products, water, steam, and chemicals
- Sanitized powder coating activity fights off bacteria and virus growth
- Remarkable graphic motifs available
- Available in galvanized steel, stainless steel, or glass

Doors

- Touchless activation option to eliminate bacteria or virus spread
- Sliding, hinged, automatic, and manual door options available

Window and Glass Elements

- Components sit flush within panel system to increase room cleaning efficiencies
- Variety of sizes and applications



Shorter
Renovation
Downtime



Hygienic
Environment



Customizable to
Your Needs



Essential Stainless Solutions

Pass Through Windows & Peg Boards

- Maintain a decontaminated environment
- Ideal for maintaining strong sterile processing department flow
- Manufactured of heavy-duty stainless steel
- Built to support constant, heavy healthcare usage

Stainless Steel Carts

- Make transportation of multiple materials easier
- Lightweight, durable, and easy to clean
- Non-porous surface
- Available in a full line of mobile opened or closed carts
- Feature fully welded construction for strength and durability

Stainless Tables

- Designed with precision engineering
- Made from high-quality materials
- Withstand the rigors of daily use
- Variety of sizes and features available

Work & Central Supply Tables

- Spacious and durable tabletop
- Ample storage with multiple shelves/drawers
- Durable, stable base for busy environments
- Ideal for organizing medical supplies, equipment, and patient records



Clean
Workspace



Customizable to
Your Needs



Skytron Assure Service Programs

A properly maintained capital asset offers years of trouble-free operation, reducing downtime and helping to control operational costs, which can improve departmental productivity. Skytron Assure provides peace of mind that your Skytron assets are performing optimally and are safe to use. Assure offers various options, including extended warranties and comprehensive service plans, all of which provide industry-leading support. Our goal at Skytron is to provide superior and affordable service to all healthcare partners and to raise the standards of our solution packages.

Advantages offered by Skytron Assure

- Maximum uptime guarantee with comprehensive plans
- Fixed or controlled operational costs
- Full life-cycle management
- Increased asset longevity
- Lowest cost of ownership over equipment life
- Support packages tailored to your needs



Maximum
Uptime
Guarantee



Prevent
Unexpected
Costs



Protect
Your
Investment

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV4 05-25
PN L2-010-202