



Spectra

Examination Light



Advanced Surgical Light Technology in an Exam Light

Built to withstand the rigors of the busiest Emergency Departments, Critical Care Units, Labor and Delivery, Clinics, Imaging, and more!



Reduced Shadows and Crisp Illumination

Spectra's seven LEDs use a Vertically Segmented Reflector Design (VSRD) to allow optimal shadow control. In addition to the brilliant spot, caregivers will enjoy the flexible positioning and drift-free stability.

Spectra Details:

- 51,000 lux at 1 meter
- Color temperature 4300°K
- Cool, color-correct light provides CRI of 96
- Optional sterilizable handle
- LED life 40,000 hours
- LED technology provides brilliant white light with no bulbs to change
- Homogeneous light spot size of 6 inches
- Durable and sealed construction for easy cleaning
- Low profile design for ease of movement and easy storage



Energy Cost
Savings



Exceptional Visibility



Durable Construction



Several Mounting Options for Different Applications


Spectra AUT1W

Wall mounted Spectra LED exam light with 75.5" reach, 180° arm rotation and 355° lighthead rotation.



Spectra AUT1S

Spectra LED exam light on mobile floor stand with 38" reach, 50° upward and 30° downward movement.



Spectra AUT1C

Ceiling mounted Spectra LED exam light with 71" reach, 360° arm rotation and infinite lighthead rotation.

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV7 03-24
PN L2-010-85