



Stainless

Scrub Sinks

Quality stainless steel scrub sinks for every clinical need

Skytron's Stainless Steel Scrub Sinks provide quality, convenience, and flexibility for every clinical need. Scrub Sinks are available with knee-kick and/or optional infrared modules for "hands-free" operation. In addition, the optional Eye Wash (EW) feature provides a central, visible location to all staff and saves valuable floor and wall space.



Single Bay with optional Infrared Module



Optional Eye Wash (-EW) safety feature



Dual Bay with optional Infrared Modules and Eye Wash



Optional Infrared Module (-IRM) conserves water usage



Triple Bay with optional Infrared Modules



The Standard Knee Kick (-MK) panel is 16" wide and 8-1/2" high and is centered to the faucet for left or right knee operation



Durable
Stainless
Construction



Fully
Configurable
Menu



Knee Kick
Inclusion

A highly-visible, deck-mounted mixing valve provides optimal temperature control. Deep-sloped sink design minimizes backsplash, while plexiglass splash dividers prevent cross-bay contamination. High quality gooseneck faucet(s) create wide, "rosette-style" spray or steady stream dependent upon user preference. Large 1.5" drain(s) promote fast and efficient water drainage for every scrub cycle.

Designed to provide the highest degree of durability and ease of use

Differentiating Features

- Solid 304 series stainless steel base
- Adjustable stream/rosette aerator — meets individual preferences for comfort and saturation of hands with water during the scrub process
- Drain in every bay reduces the potential for unwanted bacteria development (no trough design)
- One hot/cold hook-up with optional multiple thermostatic mixing valves (applicable to dual or triple bay sinks)
- Knee kick inclusion even with IR feature (dependability in event of losing power)

Optional Equipment

- Eye wash feature
- Infrared module(s)
- Laminar flow nozzles
- Additional deck mounted mixing valve
- Knee-operated soap
- Foot-operated soap
- In-Wall support system (needed when installing on plaster and drywall)
- Overhead wall shelf



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.
To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV6 03-24
PN L2-010-11