


# **An All-In-One Solution For Complete Disinfection And Compliance Management For Your Mobile Electronics**



The ElectroClave Solution is the industry's first comprehensive system for disinfection and management of your mobile electronics. Whether you are looking to manage patient facing devices or your staff's devices, Seal Shield has you covered. Configurable and adjustable to adapt to change in workflow and to the growth of your organization's devices. From 360-degree UV-C rapid disinfection to device re-imaging to complete compliance and audit reporting, the ElectroClave is the solution your organization needs.





# **Features and Benefits**

# **Disinfection**

360-degree disinfection, from rapid 60-second disinfection cycles to customizable cycles based on workflow and efficacy requirements, positive pressure dust free disinfection bays

#### **USB** Charge

10 USB ports enabling smart charging

# **Flexible Deployment**

Wall mount with locking feature, desktop mount with locking feature, interlocking and stackable to maximize space

#### **Device Agnostic**

Multi-bay configurations to allow for multiple device form factors

# **RFID Tracking**

UHF RFID tracking that integrates to ElectroClave Portal and existing RFID infrastructures within organization

# **Mobile Device Imaging/Syncing**

Over-the-wire syncing of phones or tablets through MDM

#### **Employee Badge Access**

Provide non-touch access to ElectroClave cabinets with NFC employee badges, supports both 125kHz and 13.56MHz





# **Mobile Device Management**

Gives organization visibility on mobile device to user association, mobile device physical location, department allocation and facilitates HIPAA compliance

# **Disinfection Compliance**

Gives organization ability to establish, track and audit device and user compliance

# **System Integrations**

Ability to integrate to existing systems to include RTLS, EMR/EHR, MDM

# Notifications/Alerts

Provides customizable notifications/alerts via email and SMS to simplify the management of workflows

#### **Intelligent Maintenance**

Allows an organization to ensure greatest service up-time of ElectroClave through tracking and automated scheduling of services required

# **ElectroClave Status Monitoring**

Real-time indication of ElectroClave™ operation status and disinfection status throughout the organization

#### **Disinfection Mobile App**

Gives end users awareness of disinfection status and next disinfection cycle needed

REV0 10-19 PN L2-010-181

