



Cart Wash Detergent

Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.



Cart Wash Detergent

Cleaning and Disinfection Solution

Designed exclusively for high-pressure automated cart washers, Skytron's Cart Wash Detergent is a non-foaming cleaner ideal for heavy-duty applications and deemed safe for use on aluminum surfaces. This specialized detergent is recommended for achieving spot-free drying, particularly in regions with hard water.

Recommended Uses:

- Medical carts

Advantages:

- Specially formulated for all makes of cart washing units
- Will not harm copper, stainless steel, tungsten carbide, brass or alloys
- Contains blood cutting agents and rust inhibitors
- Special surfactant blend prevents spotting
- Super concentrated to save money

Technical Data:

- pH neutral
- Non-flammable
- Contains no foaming phosphates, chlorides or hydroxides
- Bio-degradable
- Dilution – 1 oz. per 10 gallons

Packaging:

- One 5 gallon container

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 12-23
PN L2-010-241