



# Cart Wash Detergent

Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.

Scan the QR code  
to see Skytron's full  
offering of chemicals.



# Cart Wash Detergent

## Cleaning and Disinfection Solution

Designed exclusively for high-pressure automated cart washers, Skytron's Cart Wash Detergent is a non-foaming cleaner ideal for heavy-duty applications and deemed safe for use on aluminum surfaces. This specialized detergent is recommended for achieving spot-free drying, particularly in regions with hard water.

### Recommended Uses:

- Medical carts

### Advantages:

- Specially formulated for all makes of cart washing units
- Will not harm copper, stainless steel, tungsten carbide, brass, or alloys
- Contains blood cutting agents and rust inhibitors
- Special surfactant blend prevents spotting
- Super concentrated to save money

### Technical Data:

- pH neutral
- Non-flammable
- Contains no foaming phosphates, chlorides, or hydroxides
- Bio-degradable
- Dilution – 1 oz. per 10 gallons

### Packaging:

- One 5 gallon container – CSTART-8305

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV2 08-24  
PN L2-010-241