



# Cart Rinse

Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.

Scan the QR code  
to see Skytron's full  
offering of chemicals.



# Cart Rinse

## Cleaning and Disinfection Solution

Skytron's Cart Rinse is specifically formulated for use in automatic cart and spray washers. It is a highly concentrated solution tailored for application in cart and automated washer decontaminators. Additionally, Skytron's Cart Rinse effectively conditions the water, promoting virtually spot-free drying.

### Recommended Uses:

- Medical carts

### Advantages:

- Specially formulated for all makes of cart washing equipment
- Special surfactants provide supplementary cleaning
- Helps rinse away spots and film
- Creates a powerful "sheeting" action to reduce drying time

### Technical Data:

- pH 8.2 – 9.2
- Non-corrosive
- Non-flammable
- Bio-degradable
- Dilution – 1-3 oz. per 10 gallons of water

### Packaging:

- One 5 gallon container – CSTART-8105

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV2 08-24  
PN L2-010-249