



# Neutralizing Rinse

Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.



# Neutralizing Rinse

## Cleaning and Disinfection Solution

Neutralizing Rinse is a gentle acidic solution crafted to counterbalance the alkaline buildup left by detergents on various surfaces such as glassware, surgical instruments, rubber, metal, and plastic items. It effectively eliminates spots and maintains the pristine appearance of instruments and equipment. Additionally, Neutralizing Rinse guarantees thorough neutralization of any lingering alkaline detergent residue from the cleaning process in automatic washer/decontaminators.

### Recommended Uses:

- Medical/surgical instruments
- Glassware
- Rubber, metal, and plastic items
- Respiratory therapy equipment

### Packaging:

- Four 1 gallon containers – CSTART-7701
- Two 2.5 gallon containers – CSTART-7725
- One 5 gallon container – CSTART-7705

### Technical Data:

- pH 6.5 to 7.5
- Non-flammable
- Dilution – 1-3 oz. per rinse cycle

### Advantages:

- Acidic rinse promotes natural passivation process for additional instrument protection
- The anti-corrosive agents inhibit corrosion which prolongs the instruments use life
- Chelating additives counteract hard water and help protect surgical instruments by preventing spotting and staining
- Super concentrated to save money

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 01-24  
PN L2-010-251