



# 3008 Clarity

## Surgical Table

The 3008 Clarity brings Skytron quality and versatility to an imaging surgical table. With a radiolucent 83" carbon fiber table top and floating top functionality, the 3008 offers precise and infinite control of the top's movement. The Clarity's wide range of motion and radiolucency make it a versatile and adaptive table for a variety of specialties and case types.

### Imaging Access

The 3008 Clarity is equipped with an 83" carbon fiber top capable of sliding in any direction. Fully radiolucent, the carbon fiber construction allows clinicians to use the lowest doses for imaging and still maintain excellent image clarity. The 3008 boasts an impressive 65" unobstructed longitudinal imaging window.

### Floating Top

The float function is engaged with a single button, allowing the user to reposition the top in any direction with ease even with a 500lb patient load.

### Hand Control

A specially designed, comfort grip hand control is all that is needed to articulate the 3008 Clarity floating top. The float function is engaged with a single button, allowing the user to push or pull the top in any direction with ease.



# 3008 Clarity

## Surgical Table



### SPECIFICATIONS

#### Dimensions

Table Top Length	82.5" (2,100 mm)
Table Top Width	19 5/8" (500 mm)
Width w/ side rails	22" (560 mm)

#### Articulations

Height Max	41 3/4" (1065 mm)
Height Low	27 3/4" (710 mm)
Lateral Tilt (left/right)	20°
Trendelenburg	20°
Reverse Trendelenburg	20°
Slide X, Y	8" x 28"

#### Weight Capacities

Weight Cap Lift & Top Slide	500 lb.
Weight Cap Art	300 lb.

#### Pendant Control

Return to level (return to horizontal position from articulated positions)

#### Hand Console

Lateral Slide  
Longitudinal Slide

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at [www.skytron.com](http://www.skytron.com)

Single-touch pendant and free float controls

- permits smooth multi-directional top slide
- operation from any point around the table's perimeter

Full function, electro-hydraulic battery powered surgical table

- 500 lb lift and top slide capacity
- four-point, self-levelling brake system assures stability and patient safety
- automated return to level function



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REVO 12-19  
PN L2-010-185