



# 3503 EZ Slide

Surgical Table

# Simplify full body imaging and patient positioning



Featuring convenient top slide and the highest table height in our portfolio, the 3503 EZ Slide is designed with the surgical team in mind. Imaging and surgical access is easy with 21" of top slide, allowing for movement of the patient through the imaging window with just the click of a button. With a high height of 43", the 3503 EZ Slide gives comfortable working access to a broad range of specialties.

## Imaging Efficiency

With 21" of top slide, the 3503 EZ Slide makes full body imaging convenient. A simple click of a button effortlessly moves the patient toward the C-arm, eliminating the need for risky repositioning. The table also features a slim profile base that stays clear of mobile imaging equipment.

## Surgical Positioning Flexibility

The 3503 EZ Slide makes positioning and setup a breeze with a removable leg section. This can quickly integrate with accessories such as the optional 40" carbon fiber extension. Sturdy, standardized side rails accept most table accessories and preprogrammed positioning options allow for one-touch maneuvers at any point throughout the procedure.

## Improved Patient Safety

The 3503 EZ Slide takes the confusion out of weight limits. The 3503 EZ Slide can provide top slide at full weight capacity, which means the table will support articulation capacities up to 600 pounds — even while sliding the patient through the imaging window.



Comfortable  
Surgical  
Access



Powerful  
Articulation  
Capacity



One-Touch  
Imaging  
Access



## 3503 EZ Slide Series

### 21" Top Slide, General Purpose Table

Providing 21" of top slide, the 3503 EZ Slide takes the work out of imaging setup. At the click of a button, the patient is moved through the expansive imaging window.

### Differentiating Features

- 21" top slide
- 7½" slide to head and 13" slide to foot
- 700 lb lift and 600 lb articulation weight capacity
- Removable head and leg sections for quick accessory attachment and positioning
- Durable stainless steel base is easy to clean
- 25" - 43" adjustable table height
- 20° lateral tilt
- 30° Trendelenburg and 30° reverse Trendelenburg





## Convenient and Comfortable Imaging Access

The 3503 EZ Slide will offer 21" of top slide at its full articulation weight capacity. Designed to give comfortable, waist-height access to taller surgeons, the table offers our maximum table height of 43". As the table top slides toward the foot end, the head section can be conveniently detached and repositioned as a foot extension.

The 3503 EZ Slide also brings confidence to the operating room. Recognized for their reliable performance, Skytron surgical tables are consistently rated as having a lower cost of ownership.



25" Minimum  
Table Height



43" Maximum  
Table Height



21" Total Slide




40° Flex



90° Reflex



98° Leg Down



Back lit one-touch positioning functions include Trendelenburg, tilt, and more



Battery life/diagnostic indicator



Easy to use clip style power cord



Back-up controls for emergency situations



Optional 40" carbon fiber leg section for lower body imaging



Self-leveling brakes for stability on uneven floors

## One-Touch Pendant Control

- Simple one-touch operation
- Backlit control for high visibility in low lighting
- Audible alarm for kidney bridge, back up articulation and top slide obstruction
- With PSR (Patient Safe Return)
- Recessed buttons prevent accidental activation
- Durable and water resistant
- Can be used in conjunction with optional three function foot control



90° Back Up



40° Back Down



30° Trendelenburg



30° Reverse Trendelenburg



20° Tilt - Left



20° Tilt - Right

# Positioning Guide



Upper Body Imaging



Lower Body Imaging



Lap Nissen



Spinal



Ophthalmic/ENT



Nephrectomy



Neuro (Neck)



Neuro (lumbar)



Shoulder Arthroscopy



Urology



Lower Extremity



Bariatric



# Circle of Surgical Table Solutions



**3603 Ultra Slide**  
High Performance  
Removable Back  
26½" Top Slide



**6702 Hercules**  
Rotating Top &  
Removable Back



**3503 EZ Slide**  
High Performance  
Top Slide



**GS70 Salus**  
The Premier Solution  
Smart, Safe, Simple



**6302 Elite**  
Rotating Top



**3008 Clarity**  
Premier Imaging Table



**1602 Essentia**  
General Surgery Table

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at [www.skytron.com](http://www.skytron.com)



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 8-18  
PN L2-010-140