



SkyVision Ascend Eco

Video Integration

The SkyVision Ascend Eco is the next generation of medical image processing and control. Skytron has developed a powerful digital video integration solution that offers compatibility for today and tomorrow.

Users can save and recall custom presets, choose from multi-window layouts, and control image zoom from the user touch-panel interface.



Ascend Eco

What You Want to See, Where You Want to See It

The Ascend Eco provides seamless connectivity in a compact footprint. This system efficiently routes and manages all connected sources within the operating room, ensuring a streamlined surgical experience. Ascend Eco's vendor neutral design allows it to accommodate any imaging vendor or connected source with its 16 available inputs.

This cost-effective solution delivers the essential features required for seamless data routing and source management. This makes the Ascend Eco the ideal choice for facilities seeking a basic routing system without compromising performance.

Key Features

- Built-in 7" back up control touch screen
- Provides PIP, PBP, Quad, and zoom layout controls
- Live preview of selected source on touch panel control
- Can be installed remotely up to 300' from the OR
- Create, save, and edit 10 user presets
- Connect any 4K input to any 4K output
- Multi-modality routing

| | |
|----------|------------------|
| Inputs | 16 |
| Outputs | 12 |
| Supports | SD, HD, 4K Video |



Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 06-23
PN L2-010-233