



SkyVision Ascend IP Eco

Video Integration

Meet the SkyVision Ascend IP Eco, the next evolution in video routing for surgical spaces that demand simplicity, efficiency, and performance. Built on the trusted SkyVision architecture, this compact switch connects up to 10 inputs and 4 outputs effortlessly.



Smart, Space-Saving OR Control

The Ascend IP Eco is controlled through a wall-mounted 10-inch tablet with an intuitive interface that allows staff to easily manage routing and system functions. Built on SDVoE architecture, it supports seamless, vendor agnostic integration with third-party devices. Its compact, wall-mounted enclosure reduces OR clutter and space requirements, while remaining compatible with analog, HD, and 4K devices.

Space-Saving, Zero-Footprint Design

- Wall-mountable, conceal-ready, or easily integrated into existing infrastructure
- Compact, low-profile design eliminates bulky hardware
- Maximizes valuable space while maintaining a clean, uncluttered clinical environment

Intuitive Touchscreen Control

- Available in a 10-inch touchscreen for essential control of advanced layouts and routing
- Sleek, responsive interface designed specifically for surgical environments
- Pre-configured layouts and intuitive navigation enable fast, confident operation without disruption

Silent and Sustainable Operation

- No external power supplies or complex installation requirements
- Fanless, silent operation ideal for sterile environments
- Energy-efficient design supports long-term sustainability goals

