



# Sterling<sup>®</sup> C2000

Cart Washer-Disinfectors



Instrument Cart



Container Cart

## Increase Productivity with Faster Process Times

The Sterling C2000 Cart Washer streamlines and automates the decontamination process — bringing unsurpassed efficiency and flexibility to the Sterile Processing Department. The C2000 is engineered to be environmentally friendly by conserving water, energy, and chemical consumption using its built-in control functions and water recirculation system.

### Lower Cost of Ownership

- Lower water, chemical, and electric use than other washers
- Fast cycle times for increased throughput

### Efficient Work Environment


- Noise level is less than 63 dB
- Optional steam-heated drying system

### Maximize Your Space

- Four different sizes available lowers the need for custom construction
- Optional instrument cycles to augment existing washer-disinfectors



Water, Chemical,  
& Electricity  
Conservation



Low  
Noise  
Level



Fast  
Process  
Time

# Efficient Water Coverage Throughout the Chamber

- 20 wash arms with oscillating nozzles
  - **C2000** - 100 oscillating spray nozzles
  - **C2000L** - 120 oscillating spray nozzles
- Pressure controlled based on load configuration
- Wide spray pattern for complete coverage
- Spray arms and nozzles can be cleaned without tools



## The Sterling® Family of Washer-Disinfectors



The Sterling family of Washer-Disinfectors provides a full array of options for size and throughput to outfit any Sterile Processing Department with exactly what they need. Each model offers water, chemical, and electricity savings by design.

### Comply With Your Medical Device IFU's

- Edit cycle parameters to meet your needs
- Document the performance of all cycles electronically

### Boost Your Efficiency

- Best throughput per cycle in the least amount of space
- Choose the perfect sized solution for your location
- Adjustable water pressure from 9 to 14 PSI

# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Configurable Window and Door Offerings  
Built-in Cabinetry



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV3 11-24  
PN L2-010-215