












Item No.	Description	Image
S1220	Skytron Verification Test Kit - 1 vial bracket to secure in the basket, 5 ultrasonic test vials, 1 temperature strip, 1 lumen test strip holder, 5 test strips, 1 record keeping sheet and IFU	
S1225	Skytron Ultrasonic Testing Vials 25 vials and 1 record keeping sheet	
S1230	Skytron Lumen Testing Strips 20 lumen test strips, 1 record keeping sheet	
S3112	Skytron Replacement Filtering Assembly 12 in-line 70 micron filters	
S3113	Skytron Replacement Return Tubes - 6 tubes	
S3200	Skytron Male to Male Luer Slip Connectors 25 male to male luer connectors	
S3205	Latex Free Rubber Caps Small Hole for shaft diameter of 1.7mm to 5mm - 12 caps	
S3212	Latex Free Rubber Caps Medium Hole for shaft diameter of 5mm to 10mm - 12 caps	
S3215	Latex Free Rubber Caps Large Hole for shaft diameter of 10mm to 12mm - 12 caps	
S3220	Latex Free Rubber Caps Assorted Sizes (6 sm, 6 med, 6 lg)	
S3312	Large Suction Distals (1 distal)	

Item No.	Description	Image
S3315	Small Robotic Suction Distals (1 distal)	
S3316	Large Robotic Suction Distals (1 distal)	
S3320	Assorted Robotic Suction Distal Kit (3 sm distals, 1 lg distal, 3 med caps, 1 lg caps, 24 male to male connectors)	
S5000	SurgiSoak Enzyme	
CSTART-7125	Medical Enzyme Detergent Super Concentrate (2 x 2.5 gal/cs)	
CSTART-7105	Medical Enzyme Detergent Super Concentrate (1 x 5 gal/jug)	