



Cart Rinse

Clean Start Chemicals

Specifically designed for Skytron washer-disinfector systems.



Cart Rinse

Cleaning and Disinfection Solution

Skytron's Cart Rinse is specifically formulated for use in automatic cart and spray washers. It is a highly concentrated solution tailored for application in cart and automated washer decontaminators. Additionally, Skytron's Cart Rinse effectively conditions the water, promoting virtually spot-free drying.

Recommended Uses:

- Medical carts

Advantages:

- Specially formulated for all makes of cart washing equipment
- Special surfactants provide supplementary cleaning
- Helps rinse away spots and film
- Creates a powerful "sheeting" action to reduce drying time

Technical Data:

- pH 8.2 –9.2
- Non-corrosive
- Non-flammable
- Bio-degradable
- Dilution 1-3 oz. per 10 gallons of water

Packaging:

- 5 Gallons

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 12-23
PN L2-010-XXX