



Skytron

Disposable Camera Cover

Skytron's disposable camera covers are designed for use with Skytron surgical camera lights. The disposable covers are used over the camera and allow operating staff the ability to adjust the focus of the surgical light from within the sterile field. The covers come prepackaged and ready for use.

- Packaged 50 per case
- Container size - 21" W x 16" D x 19" H
- Product is sterilized using Ethylene Oxide
- Made in the USA
- Shelf life up to 5 years
- Latex free
- DEHP free
- BPA free
- No MSDS required

Skytron

Disposable Camera Cover



B1-715-65

- Simple snap-button engagement
- Single use for a Sterile Surface every procedure
- Skytron engineered to fit Stellar, Aurora LED, Aurora AUR, Aurora AUA, and Aurora AUT series light fixtures
- Clear lens area for uncompromised view and focus retention

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 5-19
PN L2-010-107