



Stainless

ES Economy Scrub Sinks

Skytron's Stainless Steel ES Scrub Sinks offer a compact design ideal for smaller spaces, delivering both quality and convenience. The ES Economy Scrub Sinks come equipped with a standard knee-kick and spout for hands-free operation, making them the perfect space-saving solution.

Skytron also offers a comprehensive line of scrub sinks (MK Series) designed to meet various functional requirements, tailored to specific project needs.

Stainless Scrub Sinks



Single Bay

ES Sink: 25"W 20-1/2"D 23-9/16"H



Dual Bay

ES Sink: 47"W 20-1/2"D 23-9/16"H



Triple Bay

ES Sink: 76"W 20-1/2"D 23-9/16"H

ES Economy Scrub Sink Features

- Single, dual, or triple stations
- Knee operated water
- Deep sloping basin
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Check mixing valve stop strainers for ease of maintenance
- Mixing valve (under sink)
- Flow control valve
- Backflow preventers
- Mounts with Z-brackets
- 304 stainless steel polished and ground to a #4 finish
- Overhead wall shelf over sink

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 12-24
PN L2-010-261