



Refer to SKYTRON Mounting Structure Pre-Installation Guide for appropriate mounting structure construction and a definition of SKYTRON's and Customer's pre-installation responsibilities.

Final model configuration requires SKYTRON technical review and approval to determine carrier configurations and equipment load capacities.




FLOOR _____ 0"

Ceiling height requirements		
	Maximum	Minimum
Carrier Fixed arm	120" [3048mm] = 72" [1829mm] top of carrier	102" [2591mm] = 70" [1778mm] top of carrier
UB Spring and elevato arm	120" [3048mm] = 78" [1981mm] bottom of UB	115" [2921mm] = 78" [1981mm] bottom of UB
Carrier Spring and elevato arm	120" [3048mm] = 71" [1803mm] center of monitor	105" [2667mm] = 70" [1778mm] top of carrier

FIXTURE WEIGHT: 518 LBS. MOMENT LOAD: 2310 FT. LBS.



Single upper arm requires bearing
H2-210-13-A



Publish Date	Revision Number	Revision History
4/20/2015	0	Initial Release
3/28/2016	1	3Q0 Hub only mount label added to page 2
7/26/2016	2	Updated dimensions