


Refer to SKYTRON Mounting Structure Pre-Installation Guide for appropriate mounting structure construction and a definition of SKYTRON's and Customer's pre-installation responsibilities.
 Final model configuration requires SKYTRON technical review and approval to determine carrier configurations and equipment load capacities.

FIXTURE WEIGHT: 628 LBS. MOMENT LOAD: 3968 FT. LBS.

Ceiling height requirements		
	Maximum	Minimum
Carrier Fixed arm	120" [3048mm] = 70" [1778mm] top of carrier	104" [2642mm] = 70" [1778mm] top of carrier
UB Spring and elevato arm	120" [3048mm] = 81" [2057mm] bottom of UB	118" [2997mm] = 78" [1981mm] bottom of UB
Carrier Spring and elevato arm	120" [3048mm] = 71" [1803mm] top of carrier	106" [2700mm] = 70" [1778mm] top of carrier



Single upper arm requires bearing
H2-210-13-A



Publish Date	Revision Number	Revision History
4/20/2015	0	Initial Release
7/26/2016	1	Dimensions corrected