



Ergostat®

Prep & Pack Tables





Personalization for Instrument Sterilization

ErgoStat® Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly, and packaging in your Department. Standing in one position for hours on end, or using a workstation designed for someone else, can lead to uncomfortable neck or back pain and unsatisfied employees. Put your faculty's well-being in mind, and choose a Prep & Pack table that easily adjusts to their needs.

Specifications	Ergostat 36 X 72 E	Ergostat 36 X 72 M	Ergostat 36 X 60 E	Ergostat 36 X 60 M	Ergostat 30 X 72 E	Ergostat 30 X 72 M	Ergostat 30 X 60 E	Ergostat 30 X 60 M
Frame Height	81"	81"	81"	81"	81"	81"	81"	81"
Width	73½"	73½"	61½"	61½"	73½"	73½"	61½"	61½"
Depth	39"	39"	39"	39"	33"	33"	33"	33"
Desk Height	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"
Height Adjustment	Power	Manual	Power	Manual	Power	Manual	Power	Manual

** All sizes available with a Corian™ or stainless steel top



Stand
or Sit



Non-Porous
Work
Surface



Easy
Accessory
Attachment

Differentiating Features

- Non-porous, solid Corian™ work surface
- 26 ½" to 41 ½" height adjustable work surface
- 43" vertical riser frame raises and lowers with the base
- Powered and manual (pin-adjustable) models available
- Heavy duty inner/outer leg extrusions for increased stability
- 750 lb dynamic and 1,000 lb static weight capacity



Optional Accessories



Akro Bin



Bin Rail,
30" or 36"




Spool Holder



Vertical CPU
Holder



Electrical Strip,
30" or 36"



Notebook Holder



Footrest,
60" or 72"



Keyboard and Mouse
Platform*



12" Post,
Two Arms



Under Worksurface
Keyboard and
Mouse Tray



Flat Screen Arm
and Mount



Magnetic Marker
Board, 18" x 30"
or 18" x 36"



Multi-function Shelf,
15" x 30" or
15" x 36"



Magnification
Light



Hanging Pedestal,
6" Box or 6" Box/Box



Lower Shelf,
18" x 60", or
18" x 72"

* Must additionally purchase 12" Post, Two Arms. PN: SKY-FSP1DMBK

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions. **To learn more, visit our website at www.skytron.com**



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV6 11-24
PN L2-010-106