



Ergostat®  
Prep & Pack Tables





## Personalization for Instrument Sterilization

ErgoStat® Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly, and packaging in your Department. Standing in one position for hours on end, or using a workstation designed for someone else, can lead to uncomfortable neck or back pain and unsatisfied employees. Put your faculty's well-being in mind, and choose a Prep & Pack table that easily adjusts to their needs.

| Specifications    | Ergostat<br>36 X 72 E | Ergostat<br>36 X 72 M | Ergostat<br>36 X 60 E | Ergostat<br>36 X 60 M | Ergostat<br>30 X 72 E | Ergostat<br>30 X 72 M | Ergostat<br>30 X 60 E | Ergostat<br>30 X 60 M |
|-------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|-----------------------|
| Frame Height      | 81"                   | 81"                   | 81"                   | 81"                   | 81"                   | 81"                   | 81"                   | 81"                   |
| Width             | 73½"                  | 73½"                  | 61½"                  | 61½"                  | 73½"                  | 73½"                  | 61½"                  | 61½"                  |
| Depth             | 39"                   | 39"                   | 39"                   | 39"                   | 33"                   | 33"                   | 33"                   | 33"                   |
| Desk Height       | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           | 26½" - 41½"           |
| Height Adjustment | Power                 | Manual                | Power                 | Manual                | Power                 | Manual                | Power                 | Manual                |

\*\* All sizes available with a Corian™ or stainless steel top



Non-Porous  
Work  
Surface



Easy  
Accessory  
Attachment

## Differentiating Features

- Non-porous, solid Corian™ work surface
- 26 ½" to 41 ½" height adjustable work surface
- 43" vertical riser frame raises and lowers with the base
- Powered and manual (pin-adjustable) models available
- Heavy duty inner/outer leg extrusions for increased stability
- 750 lb dynamic and 1,000 lb static weight capacity



## Optional Accessories



Akro Bin



Bin Rail,  
30" or 36"



Spool Holder



Vertical CPU  
Holder



Electrical Strip,  
30" or 36"




Notebook Holder



Footrest,  
60" or 72"



Flat Screen  
Keyboard Mount



Flat Screen Arm



Keyboard and  
Mouse Tray



12" Post,  
Two Arms,  
One Monitor



Magnetic Marker  
Board, 18" x 30"  
or 18" x 36"



Multi-function Shelf,  
15" x 30" or  
15" x 36"



Magnification  
Light



Hanging Pedestal,  
6" Box or 6" Box/Box



Lower Shelf,  
18" x 60", or  
18" x 72"



# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Custom Window and Door Offering  
UVC Disinfection Products



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions. **To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REVS 09-22  
PN L2-010-106