



Ergostat[®]

Prep & Pack Tables

 **SKYTRON[®]**
Reliability • Versatility • Cost Savings



Ergonomic stand or sit flexibility for your sterile processing needs

ErgoStat® Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly and packaging in your Central Sterile Supply Department. Standing in one position for hours on end, or using a workstation designed for someone else, can lead to uncomfortable neck or back pain and unsatisfied employees. Put your faculty's well-being in mind, and choose a Prep & Pack table that easily adjusts to the various physical attributes and requirements of the human body – ultimately achieving better health and employee longevity.

Specifications	Ergostat 36 X 72 E	Ergostat 36 X 72 M	Ergostat 36 X 60 E	Ergostat 36 X 60 M	Ergostat 30 X 72 E	Ergostat 30 X 72 M	Ergostat 30 X 60 E	Ergostat 30 X 60 M
Frame Height	81"	81"	81"	81"	81"	81"	81"	81"
Width	73½"	73½"	61½"	61½"	73½"	73½"	61½"	61½"
Depth	39"	39"	39"	39"	33"	33"	33"	33"
Desk Height	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"	26½" - 41½"
Height Adjustment	Power	Manual	Power	Manual	Power	Manual	Power	Manual



Non-Porous
Work
Surface



Easy
Accessory
Attachment

Differentiating Features

- Non-porous solid Corian work surface
- 26 ½" to 41 ½" height adjustable work surface
- 43" vertical riser frame raises and lowers with the base
- Powered and manual (pin-adjustable) models available
- Ergonomically correct alignment of computer monitors, lighting, bins, and pegboards
- Heavy duty inner/outer leg extrusions for increased stability
- 750 lb dynamic and 1,000 lb static weight capacity



Optional Accessories



Akro Bin



Bin Rail,
30" or 36"



Spool Holder



Vertical CPU
Holder



Electrical Strip,
30" or 36"




Notebook Holder



Footrest,
60" or 72"



Flat Screen
Keyboard Mount



Flat Screen Arm



Keyboard and
Mouse Tray



12" Post,
Two Arms,
One Monitor



Magnetic Marker
Board, 18" x 30"
or 18" x 36"



Multi-function Shelf,
15" x 30" or
15" x 36"



Magnification
Light



Hanging Pedestal,
6" Box or 6" Box/Box



Lower Shelf,
18" x 60", or
18" x 72"




Circle of Infection Prevention Solutions




Integrity Sterilizers

Increased Efficiency
Smaller Footprint



UVC Disinfection

Single Cycle, Whole-Room Disinfection
Most Comprehensive Portfolio



Processing Sinks

Ergonomically Friendly Adjustable Height
Customizable to Fit Specific Needs



ErgoStat Prep & Pack Tables

Alleviate Physical Stress
Height Adjustable Work Station



Sterling Washer-Disinfectors

Water & Chemical Conservation
Medical Device IFU Compliance



Flex Ultrasonics

Clean All Instruments in One Cycle
Thermal Disinfection

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV3 11-18
PN L2-010-106