

Flex Ultrasonic Washer-Disinfector





Clean Cannulated & Non-Cannulated in the Same Load

The Skytron Flex is an ultrasonic washer-disinfector designed specifically for orthopedic, robotic, cannulated, and non-cannulated surgical instruments.

Reduce Labor and Error

- Eliminates staff risks from manual cleaning
- Independently lab tested and validated
- Irrigates with fresh, clean water and detergent
- Meets AAMI Standards
- Utilize your current baskets and vendor trays

Ease of Use

- Two thermal disinfection cycles Cannulated and Non-Cannulated
- Reprocess mix loads in the same cycle
- Meets Intuitive Surgical Inc.'s recommended parameters for ultrasonic baths
- Built-in filtration system provides critical water final rinse

Efficient and Cost Effective

- Triple the capacity of standard ultrasonic washer-disinfectors
- Multi-step cycle cleans and disinfects in less than 35 minutes
- With its thermal disinfection rinse, Flex can be used as a backup washer-disinfector
- Reprocess 24 da Vinci instruments simultaneously





Clean Water & Detergent Flush



Uzone Final Rinse



125 Pound Weight Capacity Includes Up to Six Trays per Cycle



Differentiating Features

- Unsurpassed throughput
- Easy to operate
- No proprietary trays
- Fresh water used every cycle
- Cleans orthopedic, robotic, cannulated, and non-cannulated instruments simultaneously

Specifications	
Dimensions	52" W x 41" D x 42" H
Basin Dimensions	40" W x 29" D x 10" H
Weight Capacity	125 lbs
HF Output Power	2,000 watts
Frequency	40 kHz with a sweep of ± 2 kHz
Washer Supply	480V/3 phase 50A, 60 Hz
Watts per Gallon	48 Watts/Gallon
Final Rinse	Critical Water Final Rinse

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings Washers, Ultrasonics, and Cart Washers Prep & Pack Tables & Sterilizers



Architectural & Environmental

Modular Walls Custom Window and Door Offering UVC Disinfection Products

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



Clinical

Surgical Tables Booms & Lights Full Offering of Stainless Products



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays



5085 Corporate Exchange Blvd. SE Grand Rapids, MI 49512 1.800.SKYTRON (759-8766) email: info@skytron.com www.skytron.com



REV4 1-23 PN L2-010-150