



Flex Fit 2 & Versa Fit Stirrups

Increase your surgical team's performance with Skytron Flex Fit 2 stirrups. These industry leading lithotomy stirrups set a new standard for patient positioning in the operating room, making adjustments for gynecologic, urologic and laparoscopic procedures simple and quick.

The Flex Fit 2 allows for the control of over three axes of motion with a single trigger-grip handle. Surgeons can control abduction, adduction, and lateral leg rotation without compromising the sterile field. The trigger-grip handle is designed for better control of the patient's leg, reducing the risk of injury to both patient and surgeon.

The Versa Fit increases surgery turnaround by offering a unique-multi-boot system. Three sizes of boots range from the Standard boot up to the Bariatric HD boot that offers ample space for your largest patients.



Flex Fit 2



PN 4-090-16	Standard
Patient Weight Capacity	600 lbs
Flexible Boots	Yes
Toe Width	5.25"
Heel Width	4.25"
Calf Width	7.5"
Lithotomy Range	-25° to 75°
Adduction/Abduction Range	-9°-25°

Versa Fit



PN 4-090-17	Standard	Elite	Heavy Duty
Part Number	PN 4-090-17-1	PN 4-090-17-2	PN 4-090-17-3
Patient Weight Capacity	800	800	800
Flexible Boots	Yes	Yes	Yes
Toe Width	5.25"	5.25"	7.25"
Heel Width	4.25"	4.25"	6.25"
Calf Width	7.5"	7.5"	9.5"
Lithotomy Range	-25° to 75°	-25° to 75°	-25° to 75°
Adduction/Abduction Range	-9°-25°	-9°-25°	-9°-25°

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
 Grand Rapids, MI 49512
 1.800.SKYTRON (759-8766)
 email: info@skytron.com
www.skytron.com



REVO 11-20
 PN L2-010-208