



# Flex Fit 2 & Versa Fit Stirrups

PN 4-090-16

PN 4-090-17

**Increase your surgeon's performance** with Skytron Flex Fit 2 stirrups. These industry leading lithotomy stirrups set a new standard for patient positioning in the operating room, making adjustments for gynecologic, urologic and laparoscopic procedures simple and quick.

The Flex Fit 2 allows for the control of three degrees of motion with a single trigger-grip handle, which means surgeons can control abduction, adduction, and lateral leg rotation without compromising the sterile field. The trigger-grip handle is designed for better control of the patient's leg, reducing the risk of injury to both patient and surgeon.

## Product Highlights

- 600 lbs patient weight capacity allows you to handle more cases
- Trigger-grip handle allows for simple leg positioning
- Lithotomy and length indicators help quickly position legs
- Free-floating boots reduce the risk of superficial nerve injury
- Deluxe pads wrap around foot, ankle and calf to provide safe, comfortable and secure positioning for patients of all sizes
- Flexible boots accommodate larger patients



# Flex Fit 2 & Versa Fit Stirrups

PN 4-090-06

Built to handle a higher weight capacity while allowing simple leg positioning, the Flex Fit 2 is designed with both the patient and surgical team in mind.



## Flex Fit 2

Patient Weight Capacity	600 lbs
Flexible Boots	Yes
Toe Width	6.5"
Heel Width	4"
Calf Width	8"
Lithotomy Range	+75° to 25°
Adduction/Abduction Range	+25° to -9°

## Versa Fit

Patient Weight Capacity	800 lbs
Flexible Boots	Yes
Toe Width	6.5"
Heel Width	4"
Calf Width	8"
Lithotomy Range	+80° to -30°
Adduction/Abduction Range	+25° to -9°

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 03-24  
PN L2-010-163