



GS50 Soteria

Surgical Table

Designed for versatility, efficiency, and ease of use, the GS50 Soteria surgical table delivers the flexibility today's operating rooms demand. Advanced positioning capabilities, enhanced imaging access, intuitive controls, and broad accessory compatibility work together to support a wide range of procedures while helping clinicians deliver exceptional patient care.





Designed for Safer, Smarter Operation

Designed with safety at its core, the GS50 Soteria incorporates intuitive haptic feedback through the hand control, providing users with tactile confirmation during operation. Integrated telemetry capabilities offer valuable performance insights, while illuminated base indicators provide quick visual status updates for enhanced awareness.

Designed Around Your Needs. Expertly Designed for Longevity.



Built for Enhanced Workflow Efficiency

A flat, easy-to-clean base supports decontamination efforts and simplifies room turnover, while the waterproof pendant delivers exceptional durability in demanding surgical environments. Built for longevity and reliable performance, every detail reflects Skytron's commitment to quality, efficiency, and the real-world needs of today's operating rooms.



Exceptional Value for Everyday Surgery

Designed with intuitive controls, durable construction, and versatile positioning capabilities, the GS50 Soteria provides dependable support for a wide range of surgical procedures. By balancing cost-effectiveness with quality and safety, the GS50 Soteria offers an ideal solution for facilities seeking exceptional value without compromising patient care.

Developed through direct customer feedback, the GS50 Soteria surgical table delivers the features clinicians need most. Combining advanced functionality, enhanced safety, and exceptional durability, the GS50 Soteria builds on the trusted performance and innovation Skytron customers have come to expect from its surgical table portfolio.

Differentiating Features


- Compatible with G-Series accessories
- 8.74" (22.2cm) slide toward head and 8.74" (22.2cm) slide toward foot
- 1000lb (450kg) lift and 600lb (270kg) articulation weight capacity
- One-touch setting for beach chair and flex
- Durable stainless steel base is easy to clean
- 25° lateral tilt
- 30° Trendelenburg / Reverse Trendelenburg
- Adjustable table height
 - 25" to 42.9" (63.5cm to 109cm)
- Removable head, back and leg section
- Hand Pendant displays current position
- Reverse patient orientation mode
- Waterproof pendant with haptic feedback



Optimally Designed for Lithotomy Procedures

Purpose-built to simplify procedures and improve workflow efficiency, the table's innovative Lithotomy Mode automatically recognizes when the leg section has been removed, reducing setup steps and helping clinicians position patients more efficiently. By streamlining transitions and minimizing manual adjustments, staff can focus more on patient care and less on equipment management.





Safety at the Core of Every Procedure

- **Intuitive Haptic Feedback:** The hand control delivers tactile confirmation with every command, giving clinicians confidence that adjustments have been received and executed accurately.
- **Real-Time Telemetry:** Integrated telemetry provides continuous visibility into key table positions—including height, trend, tilt, and slide—enhancing control and situational awareness throughout procedures.
- **Instant Visual Status Indicators:** Illuminated base lights clearly display power and brake status, allowing staff to quickly verify the table is secure and ready for use.
- **Safety-Driven Design:** Together, these intelligent features improve communication between the user and the table, helping increase efficiency, and promote a safer environment for patients and surgical teams.

Intelligent Monitoring for Greater Peace of Mind




Color screen

Backlit control for high
visibility in low lighting

Recessed buttons
prevent accidental
activation with
lockout mode

Waterproof design



Skytron elevation columns
set the standard for a
solid, stable platform for
surgeons.

Updated release levers increase
safety and allow full table slide
in any position.



Easy-to-access backup
controls.

Robust stainless steel brakes
provide stability that keep your
patients and staff safe.



Positioning Guide



Upper Body Imaging
(Head Side)



Upper Body Imaging
(Leg Side)



Upper Body Imaging with
44" (112 cm) Carbon Fiber Extension



Lower Body Imaging



Lower Body Imaging with
44" (112 cm) Carbon Fiber Extension



Ophthalmic/ENT



Neuro (Neck)



Neuro (Lumbar)



Spinal



Shoulder Arthroscopy



Lap Nissen



Urology




Nephrectomy

Bariatric Positioning

1000 lb (450 kg) Lift & 600 lb (270 kg) Articulation Weight Capacity

As the demand for bariatric and higher-weight patient procedures continues to grow, having an operating table designed to safely accommodate these patients is essential. The GS50 Soteria provides the strength, stability, and positioning flexibility required to support a wider range of patient populations while helping clinicians maintain optimal access throughout the procedure.





Leading the Way in Efficiency, Positioning, and Patient Safety.

Skytron's Industry-Leading Surgical Tables

- Simplifies imaging access with intuitive controls and a low-profile base for smooth patient positioning
- Accommodates multiple surgical specialties with removable back and leg sections
- Reduces setup and turnover times with preprogrammed positioning options
- Minimizes patient handling through advanced table rotation while keeping the base securely locked
- Improves staff awareness with position timers, instability alerts, and vocalized notifications
- Enhances safety with Guardian Series smart accessory recognition to help prevent collisions

Skytron innovations build on a 50-year legacy of creating patient and staff-centered environments that prioritize safety and efficiency. We take a consultative approach, actively listening to our customers' challenges and initiatives, enabling us to align our solutions to streamline processes, reduce risk, and deliver maximum value. Our dedicated network of local specialists provides a consistent experience, supporting you from initial design through implementation and ongoing support.



5085 Corporate Exchange Blvd SE, Grand Rapids, MI 49512 | 1.800.SKYTRON (759-8766)
Email: Info@Skytron.com | www.Skytron.com



REV1 06-26
PN L2-010-283