



# Hybrid Suite

Flexible Designs





## Getting the most from your OR

Increasingly, clinicians require more versatile environments that can enable them to carry out complex minimally invasive procedures. Skytron has teamed up with different imaging vendors to create hybrid OR solutions that save clinicians critical time and flexibility to perform open or complex minimally invasive procedures in the same room.

### Versatile environments that optimize efficiency


- Easy two-touch routing with Skytron Ascend integration
- Freedom Booms offer unparalleled reach to put equipment where it needs to be

### Carry out complex, minimally invasive procedures

- Clinical staff benefit from a more tailored and intuitive environment
- Reduce costs for the care provider and decrease the amount of time the patient spends in the hospital

### Future proof solutions

- Equipment
- Plug-and-Play installation makes replacing outdated technology simple



Greater Reach  
& Vertical  
Flexibility



Minimum 360°  
Rotation



Easy to  
Move with  
Active Assist





## The perfect match of technology & experience

Developing a well-planned and flexible 'hybrid' solution for healthcare providers increases utilization of the room. Skytron's experience with interventional imaging companies provides optimal solutions for this complex environment. Clinicians benefit from Skytron's expertise by allowing us to help in the design process to confirm the correct equipment will reach where it is needed.

**Technology and experience insures that the right equipment is available for surgical and image-guided specialists including:**

- Vascular Surgeons
- Cardiothoracic Surgeons
- Interventional Cardiologists
- Neurosurgeons
- Electrophysiologists



## Expert design and support


Designing unique, complex Hybrid Rooms takes time and expertise. Skytron offers flexible solutions for designing these types of rooms, whether you are designing for today's needs or future procedures. Mounting multiple arms from a single spindle allows for greater flexibility in positioning your equipment. The Freedom Boom series provides the ability to adapt to future needs by adding or changing your equipment.

**Turn any idea into reality with advanced planning tools such as:**

- OR concept guide
- Project and planning guide
- Online portal for architects and planners
- PowerTrak 3-D room designer

# Freedom Carrier Configurations

Freedom offers many configurations to support your specific equipment and utility needs not only today but with ample space for future expansion.



**4CMC/4CMW Series**  
with front supports

- Three-point rail design
- Center-aligned shelves
- 110 lb capacity per shelf
- 37" & 49" heights
- 17", 22" or 32" shelf widths



**4CMC/4CMW Series**

- Open design
- Center-aligned shelves
- 77 lb capacity per shelf
- 37" & 49" heights
- 17", 22" or 32" shelf widths



**4VBC/4VBW Series**


- Open design
- Cantilevered shelves
- 77 lb capacity per shelf
- 37" & 49" heights
- 17", 22" or 32" shelf widths



**4UBW Series**

- Horizontal utility box
- Under-mount shelf
- 77 lb shelf capacity
- 8" height
- 22" shelf width


## Exclusive Articulating Hub is Perfect for Compact Spaces




Freedom Booms offer more versatility than the competition. With the ability to offer an articulating hub, more products can be mounted from a single arm which works perfectly in areas with tight spaces or when you need to store items out of the way when they are not in use.




# Ceiling-Mounted Imaging System Options




Compatible with Canon Alphenix™



Compatible with Philips FlexMove



Compatible with Siemens ARTIS pheno



Compatible with GE Discovery™ IGS 7

## Hygienic Work Environment

- Reduce installation times and bacterial growth by installing EASE Modular Walls into your Hybrid OR

## Complete Lighting Control


- Control the intensity and focus of your lights within the sterile field with our Aurora Four surgical light

## Easily Manage Video Information

- Improve communication and workflow within your Hybrid OR with our SkyVision Ascend
- Efficiently handle video and audio data with Skytron's CAST Media Management System

## Save Space and Reduce Clutter

- Attain convenient access to all supplies needed for patient care



## More Access, Savings, and Safety

- Reduce trip hazards by keeping hoses, cables and cords safely off the floor
- Convenient access to lights, equipment and flat-screen displays
- Lower construction costs with our single mount design

## Easy Manipulation

- Our light boom arms make adjusting radiation shields and other equipment smooth and effortless
- Increase coverage and position equipment from a single mount

## Clear and Precise Procedures

- Perform with precision using high resolution displays that offer quick video response times



## Versatile environments? Skytron has it covered.

No matter the clinical need, Skytron will work alongside planning teams to aid in the design of a basic room or the most complex. Skytron has more options and combinations for any budget or specialty.

### Find a solution for every environment or specialty

- General Operating Room
- **Hybrid Operating Room**
- C-Section Suite
- Interventional Lab
- Intensive Care Unit
- Emergency Department

# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Custom Window and Door Offering  
UVC Disinfection Products



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 09-22  
PN L2-010-210