



Lumos

Disposable Handle Cover & Camera Cover

Skytron's disposable handle and camera covers let staff adjust the Lumos lighting system without breaking the sterile field. Capacitive touch controls sterile intensity adjustment, so the surgical team can fine-tune focus and brightness mid-procedure. Covers arrive pre-packaged and ready to use, with no setup required.



Sterile Control, Zero Compromise

Skytron's disposable Handle Cover and Camera Cover let surgical staff adjust the Lumos light's position, focus, and intensity without leaving the sterile field. Capacitive touch controls respond instantly through the cover, so the team can fine-tune lighting mid-procedure without breaking scrub. Each cover ships pre-packaged and sterile, ready to use straight out of the box with no setup or reprocessing required.

Easy to Use

The intuitive design with sterile intensity control (capacitive touch) on the cover allows for effortless and non-invasive user experience for all clinical applications.

Sterile Barrier

The disposable covers enable operating room staff to adjust the surgical light focus from within the sterile field.

Cost Efficiency

With Lumos disposables single product packaging, users can easily adapt to their clinical environment without increasing cost.

- Packaged 20 per case
- Shelf life up to 5 years



B1-530-17



B1-530-18

