



# Stainless

## MK Scrub Sinks

Skytron's Stainless Steel Scrub Sinks provide quality, convenience, and flexibility for every clinical need. Scrub Sinks are available with knee-kick and/or optional infrared modules for "hands-free" operation. In addition, the optional Eye Wash (EW) feature provides a central, visible location to all staff and saves valuable floor and wall space.

Skytron also offers the ES Economy Scrub Sinks which feature a compact design with hands-free operation via a standard knee-kick and spout, combining quality, convenience, and space-saving efficiency.



**Single Bay**  
MK Sink: 32"W 27-1/2"D 37"H

**Quick Ship Option**  
SS2121-MK-SKY: 32"W 27-1/2"D 37"H  
Ready to ship sink with infrared control, eyewash, timer, dual laminar flow, blending valve, and a single bay overhead wall shelf



**Dual Bay**  
MK Sink: 64"W 27-1/2"D 37"H

**Quick Ship Option**  
SS2121-MK-SKY: 32"W 27-1/2"D 37"H  
Ready-to-ship sink with infrared control, eyewash, timer, dual laminar flow, blending valve, thermostatic mixing valve, and a dual bay overhead wall shelf



**Triple Bay**  
MK Sink: 96"W 27-1/2"D 37"H



# Stainless Scrub Sinks

## MK Scrub Sink Features

### Standard

- Single, dual, or triple stations
- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve (one per station)
- Knee operated water
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Splash screen dividers with dual and triple station sinks
- Main support and Z-brackets
- 304 stainless steel polished and ground to a #4 finish

### Optional

- Infrared operation
- Knee operated soap
- Foot operated soap
- Digital scrub timer
- Eye wash
- Blending valve “tepid water” for eye wash
- Additional thermostatic mixing valve
- Laminar faucet diffuser
- 220/240V power supply for infrared control (International only)
- Pedestal base
- In-wall carrier
- Overhead wall shelf over sink

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 12-24  
PN L2-010-262