



Sterling P1000 & M1000

Washer-Disinfectors



Worry-free compliance while saving time, money, and space

The Sterling M1000 & P1000 washer-disinfectors streamline and automate the decontamination process – bringing unsurpassed efficiency and flexibility to the Sterile Processing Department.

Lower Cost of Ownership

- 64% less water and chemical use than other washers
- Lowers utility and sewage bills

Comply With Your Medical Device IFU's

- Allows for 40 fully customizable cycles
- Intuitive touchscreen controls allow for editing parameters for pre-wash, wash, rinse or dry phases

Maximize Your Efficiency


- 50% more trays in 80 inches of space width
- Up to 12 DIN trays per cycle



Water
Conservation



Customizable
Cycles



Conserve
Space



Powered Door



Manual Door

Automatically select and start the next cycle when the previous cycle has finished with optional RFID equipped manifold racks

Accommodate instrument trays of any size M1000's modular manifold racks are easily configured

Guarantee consistent performance the on screen checklist for daily maintenance prompts automatic quality assurance

Protects against water quality changes utilizing an optional conductivity sensor

Configured for Versatility


Instruments



Utensils



MIS Specific



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Custom Window and Door Offering
UVC Disinfection Products



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 10-22
PN L2-010-175