

Sterling P1000 & P1800

Washer-Disinfector Accessory Guide

Setting the standard for sterile processing

The Sterling P1000 and P1800 Washer-Disinfector series streamlines and automates the decontamination process – bringing unsurpassed efficiency and flexibility to the Sterile Processing Department.

Table of Contents	1
Wash Racks	1
Transport Carts	5
Rack Inserts and Support	7
Injectors and Nozzles	9

Wash Racks

Flexible Rack

Three Levels

- Wash in up to three levels
- Includes grids and detachable washing arms
- Maximum height of 5.1 inches (13 cm) on the top level and 7.8 inches (20 cm) on each lower level
- 2 removable spray arms and shelves

Item #1000-FLX3 Holds up to 6 DIN trays

Item #1800-FLX3 Holds up to 9 DIN trays



Flexible Rack

Four Levels

- Wash in up to four levels
- Includes grids and detachable washing arms
- Maximum height of 4.6 inches (11.8 cm) on the top level and 4.6 inches (11.8 cm) on each lower level
- 3 removable spray arms and shelves

Item #1000-FLX4 Holds up to 8 DIN trays

Item #1800-FLX4 Holds up to 12 DIN trays



Flexible Rack

Five Levels

- Wash in up to five levels
- Includes grids and detachable washing arms
- Maximum height of 3.2 inches (8.2 cm) between levels
- 4 removable spray arms and shelves

Item #1000-FLX5 Holds up to 10 DIN trays

Item #1800-FLX5 Holds up to 15 DIN trays



Wash Racks

Fixed Rack

Four Levels

- The loading space in each level is up to 4.5 inches (11.5 cm)

Item #1000-FIX4 Holds up to 8 DIN trays

Item #1800-FIX4 Holds up to 12 DIN trays



Fixed Rack

Five Levels

- The loading space in each level is up to 3.15 inches (8 cm)

Item #1000-FIX5 Holds up to 10 DIN trays

Item #1800-FIX5 Holds up to 15 DIN trays



Fixed Rack

Six Levels

- Maximum height of 2 inches (5cm) at top level, max 2.4 inches (6cm) all lower levels

Item #1000-FIX6 Holds up to 12 DIN trays

Item #1800-FIX6 Holds up to 18 DIN trays



Flexible Rack Plus

5 Levels

Designed to facilitate the efficient loading and cleaning of up to 10 DIN baskets, this system features four M8x1 outlets for connecting luer or tube fittings to irrigate up to 4 instruments. Four removable spray arms and shelves for added convenience.

Loading Heights:

- All levels: 80mm
- Gain from removing spray arm and shelf: 60mm

Item #1000-FLX5-PLUS Holds up to 10 DIN trays

For P1000 or M1000 only.



Wash Racks

Flexible Rack Plus

4 Levels

Designed to facilitate the efficient loading and cleaning of up to 8 DIN baskets, this system features four M8x1 outlets for connecting luer or tube fittings to irrigate up to 4 instruments. Three removable spray arms and shelves for added convenience.

Loading Heights:

- All levels: 115mm
- Gain from removing spray arm and shelf: 60mm

Item #1000-FLX4-PLUS Holds up to 8 DIN trays

For P1000 or M1000 only.



Flexible Rack Plus

5 Levels

Designed to facilitate the efficient loading and cleaning of up to 15 DIN baskets, this system features six M8x1 outlets for connecting luer or tube fittings to irrigate up to 6 instruments. Four removable spray arms and shelves for added convenience.

Loading Heights:

- All levels: 80mm
- Gain from removing spray arm and shelf: 60mm

Item #1800-FLX5-PLUS Holds up to 15 DIN trays

For P1800 only.



Flexible Rack Plus

4 Levels

Designed to facilitate the efficient loading and cleaning of DIN baskets. The system features six M8x1 outlets for connecting luer or tube fittings to irrigate up to 6 instruments. Three removable spray arms and shelves for added convenience.

Loading Heights:

- All levels: 115mm
- Gain from removing spray arm and shelf: 60mm

Item #1800-FLX4-PLUS Holds up to 12 DIN trays

For P1800 only.



Wash Racks

Small Container Rack

The vertical container load rack is designed to wash a number of containers within the wash rack footprint.

For small containers sized 600 x 300 x 150 mm

Item #1000-CNTR-SM Holds up to 4 containers

Item #1800-CNTR-SM Holds up to 6 containers



Medium Container Rack

The vertical container load rack is designed to wash a number of containers within the wash rack footprint.

For medium containers sized 600 x 300 x 300 mm

Item #1000-CNTR-MED Holds up to 2 containers

Item #1800-CNTR-MED Holds up to 4 containers



MIC Rack

Basic Nozzles

MIC rack complete with grommet and basic nozzles. Please see Caps For Basic Nozzles on Page 11 for more information.

Item #1000-MIC-WN Comes with 40 basic nozzles

Item #1800-MIC-WN Comes with 56 basic nozzles



Wash Racks & Transport Carts

MIC Rack

Without Nozzles

MIC rack without grommet and nozzles. Please see Caps For Basic Nozzles on Page 11 to order.

Item #1000-MIC-WON

Item #1800-MIC-WON



da Vinci Rack

da Vinci rack without holders Si, Xi.

The holder capacity of this rack varies. The P1800 rack example shown is fitted with 5 Si holders, 6 Xi holders, 1 Xi Optics holder, and an RFID reader.

Item #1000-DVXISI Holds up to 10 da Vinci instrument holders or up to 6 da Vinci instrument holders and 2 da Vinci Stapler holders, and a shelf to hold 2 DIN trays

Item #1800-DVXISI Holds up to 14 da Vinci instrument holders or up to 10 da Vinci instrument holders and 2 da Vinci Stapler holders, and a shelf to hold 3 DIN trays



Manual Transport Cart

Manual loading cart transports Sterling P1000 wash racks in all configurations. Easily maneuverable around all environments.

- Fits all Skytron Sterling P1000 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1000 washer

Item #1000-MC



Transport Carts

Manual Transport Cart

Manual loading cart transports Sterling P1800 wash racks in all configurations. Easily maneuverable around all environments.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer

Item #1800-MC



Automatic Loading Cart

The Automatic loading system for the Sterling P1800 can be connected to the clean and unclean side of the machine.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer
- Minimize heavy handling of instruments

Item #1800-ALT



Automatic Unloading Cart

The Automatic loading system for the Sterling P1800 can be connected to the clean and unclean side of the machine.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer
- Minimize heavy handling of instruments

Item #1800-AULT



Rack Inserts and Support

Shelf Level Rack

Washing Stainless Bowls

Shelf level for Sterling P1000 and P1800 wash racks for washing and disinfecting kidney bowls.

- Capable of holding up to six bowls
- Can hold bowls up to a max diameter of 13.8 inches (35 cm)

Item #1000-FLX-KB



Shelf Level Rack

Washing Stainless Bowls

Fits in all flexible racks.

- Maximum diameter of 14.2 inches (36 cm) per bowl
- Capacity of 4 wash bowls
- 3 inserts per level

Item #1000-FLX-BWL Holds up to two bowl holders per shelf



da Vinci Xi Holder

da Vinci Xi Staple Holder

Fits on M/P1000 & P1800 da Vinci racks.

For optimal loading and cleaning of Xi EndoWrist Stapler instruments with two flush ports that connect to the back of the instrument housing.

- M/P1000 Able to have up to 2 da Vinci Stapler holders and 6 da Vinci instrument holders on the da Vinci rack
- P1800 Able to have up to 2 da Vinci Stapler holders and 10 da Vinci instrument holders. on the da Vinci rack

Item #1000-XIST



da Vinci Xi Holder

da Vinci Xi holder 8 mm

Fits on M/P1000 & P1800 da Vinci racks. For optimal loading and cleaning of Xi EndoWrist 8mm instruments with two flush ports that connect to the back of the instrument housing.

- M/P1000 rack holds up to 10 da Vinci instrument holders or up to 6 da Vinci instrument holders and 2 da Vinci Stapler holders
- P1800 rack holds up to 14 da Vinci instrument holders or up to 10 da Vinci instrument holders and 2 da Vinci Stapler holders

Item #1000-XI8



Rack Inserts and Support

da Vinci Si Holder

da Vinci Single Site Holder

Fits on M/P1000 & P1800 da Vinci racks. For optimal loading and cleaning of Si EndoWrist instruments with two flush ports that connect to the back of the instrument housing.

- M/P1000 Able to have up to 2 da Vinci Stapler holders on the da Vinci rack and 6 da Vinci instrument holders on the da Vinci rack
- P1800 Able to have up to 2 da Vinci Stapler holders and 10 da Vinci instrument holders on the da Vinci rack

Item #1000-SS



da Vinci Si Holder

da Vinci Si holder 8/5 mm

Fits on M/P1000 & P1800 da Vinci racks. 10 per rack.

- M/P1000 Able to have up to 2 da Vinci Stapler holders and 6 da Vinci instrument holders on the da Vinci rack
- P1800 Able to have up to 2 da Vinci Stapler holders and 10 da Vinci instrument holders on the da Vinci rack

Item #1000-SI85



da Vinci Si Holder

da Vinci Single Site Stapler Holder

Fits on M/P1000 & P1800 da Vinci racks. For optimal loading and cleaning of Si EndoWrist Single-Site instruments with primary flush port on the instrument shaft and secondary flush port on the back of the instrument housing.

- M/P1000 holds up to 2 da Vinci Stapler holders or up to 6 da Vinci instrument holders
- P1800 holds up to 2 da Vinci Stapler holders or up to 10 da Vinci instrument holders

Item #1000-SIST




Rack Inserts and Support

Fine Mesh Grate For Racks

P1000 Item #2908469

P1800 Item #2918469



Cleaning Kit

Kit for cleaning power tools.

Item #2908314



Injectors and Nozzles


Fine Filter Module

For Dental and Ophthal Units

Fine filtered module for 5 level racks for the cleaning of delicate instruments like Ophthalmic and Dental.

P1000 Item #2908539


P1800 Item # 2918539



Dental Unit 016/020

- Four - 016
- Four - 020
- Three units per module

Item #2908508



Injectors and Nozzles

Dental Unit 016

- Eight piece 016 adapter
- Three units per module

Item #2908432



Distributor Unit

Adapters

- Eight piece 020 adapters
- Two unit per module.

Item #2908506



Injector Spikes 04

Anesthesia Rack

Injector spikes for tubes 04 for Anesthesia rack.

Item #2908255



Dental Unit 020

- Four - 016
- Four - 020
- Three - units per module

Item #2908426



Injectors and Nozzles

Nozzle 04

MIC Rack

Injector nozzle screw 04 for old MIC rack and module with screw.

Item #2908364



Injector Nozzle 02.6

MIC Rack

Injector nozzle 02.6 (grommet) length 50 mm. for MIC rack and modules.

Item #2908441



Injector Nozzle Closed

MIC Rack

Injector nozzle 04 closed top (grommet) length 120mm. for MIC rack and modules.

Item #2908426



Injector Nozzle Open

MIC Rack

Injector nozzle 04 open top (grommet) length 120mm. for MIC rack and modules.

Item #2908428



Injectors and Nozzles

Basic Nozzle

MIC Rack

Basic nozzle screw for tubular instruments for MIC rack and module with screw.

Dimensions 10 x 6 x 2 mm.

Item #2900469



Basic Nozzle - Long

MIC Rack

Basic nozzle, long, (grommet) for tubular instruments.

- Maximum Ø13 mm.
- Sleeve length 55 mm.

Item #2900807



Basic Nozzle - Short

MIC Rack

Basic nozzle, short, (grommet) for tubular instruments.

- Maximum Ø13 mm.
- Sleeve length 27 mm.

Item #2900808



Ophthal Nozzle

Ophthal nozzle fine mesh filter

Item #2908373



Injectors and Nozzles

Caps For Basic Nozzles


1: Caps for basic nozzle 02 - Item #2900804

2: Caps for basic nozzle 03 - Item #2900805

3: Caps for basic nozzle 04 - Item #2900806

4: Caps for basic nozzle 03-08 - Item #2900802

5: Caps for basic nozzle 06-011 - Item #2900803



Universal Nozzle

MIC Rack

Nozzle universal screw for tubular instruments for MIC rack and module with screw.


Item #2908184



Dental Nozzle

Dental nozzle 016


Item #2908303



Dental Nozzle 020

Dental nozzle 020

Item #2908304



Nozzle 04 Closed

MIC Rack

Injector nozzle screw 04 closed for MIC rack and module with screw.

Item #2908283



Injectors and Nozzles

Nozzle 02.6

MIC Rack

Injector nozzle screw 02.6 for MIC rack and module with screw.

Item #2908439



Ophthal Unit

Three units per module

Item #2908519



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV3 07-24
PN L2-010-198