

# Sterling P1800 Washer-Disinfector





## Process your medical devices the most efficient, and compliant way possible

The Sterling P1800 is designed to set a new standard for Sterile Processing Departments. With a large capacity, small footprint, fast cycles, low water and electrical consumption and most important, disinfected devices on a consistent basis.

## Lower Cost of Ownership

- Less water and chemical use than other washers
- Save 50% in chemical and sewage costs

## Comply With Your Medical Device IFU's

- Edit cycle parameters to meet your needs
- Document the performance of all cycles electronically

## Maximize Your Efficiency

- Best throughput per cycle in the least amount of space
- Up to 18 DIN trays per cycle

Water Conservation









**Automate your SPD** with RFID equipped manifold racks that select and start your next cycle.

**Accommodate instrument trays of any size** P1800's modular manifold racks are easily configured.

**Maintain consistent performance** with a daily on-screen user checklist and remote service connectivity.

**Protects against water quality changes** utilizing an optional conductivity sensor.

**Elevate your disinfection process** with optional preheat boosters to shorten cycles and recirculation tanks to save critical water.

## Configured for Versatility

Instruments



Utensils



MIS Specific



## Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings Washers, Ultrasonics, and Cart Washers Prep & Pack Tables & Sterilizers



#### Clinical

Surgical Tables Booms & Lights Full Offering of Stainless Products



#### Architectural & Environmental

Modular Walls Custom Window and Door Offering UVC Disinfection Products

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



### Clinical Business Intelligence

Video Integration Asset Tracking & Supply Chain Solutions Surgical Displays



5085 Corporate Exchange Blvd. SE Grand Rapids, MI 49512 1.800.SKYTRON (759-8766) email: info@skytron.com www.skytron.com



REV1 09-22 PN L2-010-196