



Pediatric

Table Extension

Closing the gap between a surgeon and their littlest patients

Skytron's pediatric extension narrows the working area and improves surgeon access for smaller patients – transforming a general purpose surgical table to a pediatric specialty table.

Quick Installation

- The lightweight table extension easily attaches to the leg pins maintaining swift turnover times between cases

360° Imaging Access

- The carbon fiber construction, detachable side rails and open base combine to provide unsurpassed 360° imaging capability

Simple Accessory Attachment

- Side rails can be adjusted to any position for convenient attachment of accessories or storage for the table's pendant control



Pediatric Table Extension

Compatible with the Skytron 6300, 6500, 6600, 6700, 3500 and 3600 table series






Part Number Description

Part Number	Description
4-030-42	Pediatric Table Extension
4-030-42-2WS	2" Welded Seam Pad
4-030-42-3WS	3" Welded Seam Pad
4-030-42-4WS	4" Welded Seam Pad
4-030-42-5WS	5" Welded Seam Pad
4-030-42-9	Clamp on Rail, Pair

Key Features

- Complete 360° imaging access
- Narrow 15" wide working area
- Adjustable and detachable side rails (optional)
- Lightweight carbon fiber construction (24.6 lb)
- Attaches to the leg section of most Skytron tables
- Available with a 2", 3", 4" or 5" pad (sold separately)
- 100 lb maximum weight capacity



The literature has been developed to provide healthcare professionals with general product information. Please note: not all products are available in all countries. Product labeling and associated claims may differ based on regulatory requirements. For product availability or any further information, please contact your local Skytron representative.

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 03-24
PN L2-010-132