



Stainless

Scrub Sinks

Quality stainless steel scrub sinks for every clinical need

Skytron's Stainless Steel Scrub Sinks provide quality, convenience and flexibility for every clinical need. Scrub Sinks are available with knee-kick and/or optional infrared module for "hands-free" operation. In addition, the optional Eye Wash (EW) feature provides a central, visible location to all staff, and saves valuable floor and wall space.



Single Bay with optional Infrared Module



Optional Eye Wash (-EW) safety feature.



Double Bay with optional Infrared Modules and Eye Wash




Optional Infrared Module (-IRM) conserves water usage.



Triple Bay with optional Infrared Modules



The Standard Knee Kick (-MK) panel is 16" wide and 8-1/2" high and is centered to the faucet for left or right knee operation.



Durable
Stainless
Construction



Fully
Customizable
Menu



Knee Kick
Inclusion

A highly visible, deck mounted mixing valve provides optimal temperature control. Deep sloped sink design minimizes backsplash, while plexiglass splash dividers prevent cross-bay contamination. High quality gooseneck faucet(s) create wide, "rosette-style" spray or steady stream dependent upon user preference. Large 1-1/2" drain(s) promote fast and efficient water drainage for every scrub cycle.

Designed to provide the highest degree of durability and ease of use

Differentiating Features

- Solid 300 series Stainless Steel base
- Adjustable Stream/Rosette Aerator – meets individual preferences for comfort and saturation of hands with water during the scrub process
- Drain in every bay reduces the potential for unwanted bacteria development (no trough design)
- One hot/cold hook-up with optional multiple thermostatic mixing valves (applicable to dual or triple bay sinks)
- Knee kick inclusion even with IR feature (dependability in event of losing power)
- Original provider of eye wash feature in scrub sinks using tepid water


Optional Equipment

- Eye Wash Feature
- Infrared Module(s)
- Wall Support System
- Back Shelf
- Laminar Flow Nozzles
- Additional deck mounted mixing valve





Circle of Clinical Solutions



Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
 Grand Rapids, MI 49512
 1.800.SKYTRON (759-8766)
 email: info@skytron.com
www.skytron.com



REV3 11-18
 PN L2-010-11