



## UV Smart

### Countertop Rapid UVC Disinfectant

Achieve ultimate UVC disinfection of your personal belongings in as little as 25 seconds. The UV Smart is a simple-to-use solution to disinfect your high-touch items like cellphones, tablets, laptops, keyboards, and so much more.

 **SKYTRON**<sup>®</sup>  
Reliability • Versatility • Cost Savings

# UV Smart

## Disinfect Items Exclusive to Your Needs

### Keep Germs and Viruses Away from Patients

360-degree disinfection on high-touch items like mobile devices, remotes, keyboards, and more.

### Simple Addition to any Workflow

User-friendly design and portable making it easy to integrate into your existing space.

### Peace of Mind

Ensures a consistent and powerful disinfection in 25 seconds.



### Healthcare

- Tablets
- Phones
- Name badges
- Glasses

### Correctional Facilities

- Badges
- Radio
- Handcuffs
- Firearms

### Office Spaces

- Laptops
- Keyboards & Mice
- Cellphones
- Office Supplies

### Specifications

**Exterior Dimensions:** 21" L x 16" W x 10" H (535 L x 400 W x 260 H mm)

**Interior Dimensions:** 17" L x 10" W x 6" H (420 L x 265 W 160 H mm)

**Disinfection Time:** 25 Seconds

**Reduction:** Log-4 (99.99%)

**Disinfectant Treatment Against:** Microorganisms such as bacteria, viruses, spores, yeasts, fungi, and mycobacteria

Skytron is the Healthcare Efficiency Specialist, providing full-room solutions of capital equipment, architectural and real-time information systems for Medical, Surgical, Sterile Processing, and Infection Prevention. Our solutions enhance the utilization of people, facilities, and capital because they are designed with the user in mind and have a low, long-term cost of ownership.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REVO 03-21  
PN L2-010-216