



## UV Smart ECD

### Countertop Rapid UVC Disinfecter

Achieve ultimate UVC disinfection of your personal belongings in as little as 25 seconds. The UV Smart is a simple-to-use solution to disinfect your high-touch items like cellphones, tablets, laptops, keyboards, and so much more.



# UV Smart ECD

## Disinfect Items Exclusive to Your Needs

### Keep Germs and Viruses Away from Patients

360-degree disinfection on high-touch items like mobile devices, remotes, keyboards, and more.

### Simple Addition to any Workflow

User-friendly design and portable making it easy to integrate into your existing space.

### Peace of Mind

Ensures a consistent and powerful disinfection in 25 seconds.



### Healthcare

- Tablets
- Phones
- Name badges
- Glasses

### Correctional Facilities

- Badges
- Radio
- Handcuffs
- Firearms

### Office Spaces

- Laptops
- Keyboards & Mice
- Cellphones
- Office Supplies

### Specifications

**Exterior Dimensions:** 21" L x 16" W x 10" H (535 L x 400 W x 260 H mm)

**Interior Dimensions:** 17" L x 10" W x 6" H (420 L x 265 W 160 H mm)

**Disinfection Time:** 25 Seconds

**Reduction:** Log-4 (99.99%)

**Disinfectant Treatment Against:** Microorganisms such as bacteria, viruses, and fungi

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV2 09-22  
PN L2-010-216